

Product Datasheet

POU3F1 antibody (orb19423)

Description

Goat polyclonal antibody to POU3F1

Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Target	OCT6 / POU3F1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Preliminary experiments in lysates of Human Brain (Cerebellum, Hippocampus, Frontal Cortex), Mouse Brain, Rat Brain and Rat Spinal Cord gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	GHPPMDDVYAPGE
Clonality	Polyclonal
MW	45.5
Entrez	5453
Dilution Range	ELISA: 1:2000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.