


Product Datasheet

SMN1 / SMN2 antibody (orb19429)

Description	Goat polyclonal antibody to SMN1 / SMN2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	SMN1 / SMN2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 26kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 28.5kDa according to NP_075012.1 and NP_075014.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	DESENSRSPGNKSDN
Clonality	Polyclonal
MW	28.5; 31.8 / 30.5; 28.5; 27.1
Entrez	6606, 6607
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb19429
(1µg/ml)
staining of
Human
Cere...