


Product Datasheet

GPX3 antibody (orb19430)


Description	Goat polyclonal antibody to GPX3
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	GPX3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:1000.WB: Approx 26kDa band observed in Mouse Kidney and Liver lysates (calculated MW of 27.7kDa according to Mouse NP_001077398.1 [isoform1], 25.4kDa according to Mouse NP_032187.2 [isoform 2] and 25.5kDa according to Human NP_002075.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	NQFGKQEPGENSE
Clonality	Polyclonal
MW	25.5
Entrez	2878
Dilution Range	ELISA: 1:1000, WB: 1-3 µg/ml, IHC-P: 20 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Heart lys...



Immunofluorescence analysis of Porcine K...



Immunohistochemical staining of Human Pl...