

Product Datasheet

TRAF2 antibody (orb19446)

Description

Goat polyclonal antibody to TRAF2

Species/Host

Goat

Reactivity

Bovine, Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, IF, IHC, WB

Applications

Target

TRAF2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 50kDa band observed in Human Ovary lysates (calculated MW of 55.9kDa according to NP_066961.2). In transfected HEK293 transiently expressing Human TRAF2 a band of approx. 60kDa is observed. This band is not observed in the non-transfected HEK293. The endogenous 48kDa protein is observed in both. Recommended concentration: 0.1-0.3 µg/ml.

Protein Sequence

KMEAKNSYVRDD

Clonality

Polyclonal

MW

55.9

Entrez

7186

Dilution Range

ELISA: 1:128000, WB: 0.1-0.5 µg/ml, IHC-P: 3.75 µg/ml


Expiration Date

12 months from date of receipt.


250kDa
150kDa
100kDa

50kDa
37kDa
25kDa
20kDa
15kDa

Western blot analysis
of human ovary lys...



Immunohistochemical
staining of Human
Ki...



Immunohistochemical
staining of Human
Pl...