



Product Datasheet

Cyld antibody (orb19453)

Description	Goat polyclonal antibody to Cyld
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Cyld (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 125kDa band observed in mouse and rat brain lysates (calculated MW of 107kDa according to mouse NP_775545.1). Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	DSQQSKSKNPWYIDE
Clonality	Polyclonal
MW	107
Dilution Range	ELISA: 1:32000, WB: 0.5-1.5 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb19453 (0.5ug/ml)
staining of rat brain...



Immunohistochemical
staining of Mouse
Br...