

Product Datasheet

p21 Antibody (orb1945303)

Description	p21, also known as p21 ^{Cip1} , p21 ^{Waf1} , Cyclin-Dependent Kinase Inhibitor 1, or CDK-Interacting Pr...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Preservatives	Tris Buffer, pH 7.3 - 7.7, with 1% rAlbumin and 0.1% Sodium Azide
Storage	Store at 2-8°C. To ensure stability, immediately replace vial back in the refrigerator after each use
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.