

Product Datasheet

BTBD14B antibody (orb19457)

Description	Goat polyclonal antibody to BTBD14B
Species/Host	Goat
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	NAC1 / BTBD14B
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx 55kDa band observed in Human Bone Marrow lysates (calculated MW of 57.3kDa according to NP_443108.1). Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	KTEQQESDSVQC
Clonality	Polyclonal
MW	57.3
Entrez	112939
Dilution Range	ELISA: 1:2000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



orb19457
(0.5ug/ml)
staining of
Human Bo...