



Product Datasheet

KLK2 antibody (orb19459)

Description	Goat polyclonal antibody to KLK2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	KLK2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 22kDa band observed in Human Prostate lysates (calculated MW of 24.7kDa according to NP_001002231.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	KHQSLRPDEDSSHD
Clonality	Polyclonal
MW	28.7; 24.7;
Entrez	3817
Dilution Range	ELISA: 1:64000, WB: 0.5-2 µg/ml, IHC-P: 2.5-8µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Prostate tissue showing KLK2 expression.



Immunohistochemical staining of Human Prostate tissue showing KLK2 expression.



Immunohistochemical staining of Human Prostate tissue showing KLK2 expression.