



Product Datasheet

SARS-CoV-2 BA.5 Omicron Variant Recombinant Spike Trimer His Tag (orb1945941)

Description	SARS-CoV-2 BA.5 Omicron Variant Recombinant Spike Trimer His Tag
Reactivity	Other
Tested Applications	ELISA
Preservatives	1X PBS, 0.025% NaN ₃
Concentration	0.5 mg/ml
Storage	2-8°C
Tag	His tag
Note	For research use only
Isotype	SARS-CoV-2
Dilution Range	Control antigen for SARS-CoV-2 (COVID-19) diagnostic assays. Contains C-terminal His tag used for purification.
Expiration Date	6 months from date of receipt.



Concentration-response curves for bind...



SARS-CoV-2 Spike Trimer, His Tag (BA.5/O...