



Product Datasheet

SARS-CoV-2 Mu Variant Recombinant Spike RBD His Tag (orb1945958)

Description	SARS-CoV-2 Mu Variant Recombinant Spike RBD His Tag
Reactivity	Other
Tested Applications	ELISA
Preservatives	1X PBS, 0.025% NaN ₃
Concentration	0.5 mg/ml
Storage	2-8°C
Tag	His tag
Note	For research use only
Isotype	SARS-CoV-2
Dilution Range	Control antigen for SARS-CoV-2 (COVID-19) diagnostic assays. Contains C-terminal His tag used for purification.
Expiration Date	6 months from date of receipt.



SARS-CoV-2
Mu Spike
Trimer
Variant His
T...



SARS-CoV-2
Mu B.1.621
Variant
Recombinan...