

Product Datasheet

SARS-CoV2 G10xA5 human neutralizing mAb Antibody (orb1945963)

Description	SARS-CoV2 G10xA5 human neutralizing mAb Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FA
Immunogen	A recombinant SARS-CoV-2 spike protein containing RBD.
Preservatives	1X PBS, 0.2% rAlbumin
Concentration	0.5 mg/ml
Storage	-20°C
Note	For research use only
Isotype	Human IgG1k
Clonality	Monoclonal
Dilution Range	Experimental conditions may be optimized.
Expiration Date	12 months from date of receipt.