

Product Datasheet

SARS-CoV2 LxC1-G10 human neutralizing mAb Antibody (orb1945968)

Description

SARS-CoV2 LxC1-G10 human neutralizing

Reactivity

Other

Conjugation

Unconjugated

Tested Applications

ELISA, FA

Immunogen

A recombinant SARS-CoV-2 spike protein containing RBD.

Preservatives

1X PBS

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Human IgG1k

Clonality


Monoclonal

Dilution Range


Experimental conditions may be optimized.

Expiration Date


12 months from date of receipt.



Spike protein trimers were coated at 2 µg/ml



Spike protein trimers were coated at 2 µg/ml



Spike protein trimers were coated at 2 µg/ml