

Product Datasheet

SARS-CoV2 LxC1-E7 human mAb Antibody (orb1945970)

Description

SARS-CoV2 LxC1-E7 human mAb Antibody

Reactivity

Other

Conjugation

Unconjugated

Tested Applications

ELISA, FA

Immunogen

A recombinant SARS-CoV-2 spike protein containing RBD.

Preservatives

1X PBS

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality


Monoclonal

Dilution Range


Experimental conditions may be optimized.

Expiration Date

12 months from date of receipt.



Spike protein trimers were coated at 2 µg/ml



Spike protein trimers were coated at 2 µg/ml



SARS-CoV2 LxC1-E7 human neutralizing mAb