


Product Datasheet


SARS-CoV-2 spike protein (RBD) neutralizing (C-A11) rabbit mAb Antibody (orb1945971)

Description


SARS-CoV-2 spike protein (RBD)



Flow
cytometry
neutralizing
assay: C-
A11...



Affinity
assessment
of SARS-
CoV-2 C-
A11 ...



Flow
cytometry
neutralizing
assay of C-
A1...

Reactivity

Other

Conjugation

Unconjugated

Tested

FA

Applications

Immunogen

recombinant protein (RBD)

Preservatives

1X PBS

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality

Monoclonal

Uniprot ID

P0DTC2

Dilution Range

This antibody recognizes the SARS-CoV-2 Spike Protein RBD domain and inhibits the interaction between SARS-CoV-2 RBD and ACE2. Can be paired with other antibody for detection in sandwich ELISA format.

Expiration Date

12 months from date of receipt.