

Product Datasheet

Anti-PLXDC2 (4G3) rabbit mAb Biotin Conjugate (orb1945972)

Description

Anti-PLXDC2 (4G3) rabbit mAb Biotin Conjugate

ELISA

Reactivity

Human

Conjugation

Biotin

Tested

WB

Applications

Immunogen

HEK293 cells expressing human PLXDC2

Preservatives

1X PBS, 0.05% (w/v) Proclin300

Concentration

0.5 mg/ml

Storage

2-8°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality

Monoclonal

Uniprot ID


Q6UX71

Dilution Range

ELISA 1:5000-10000 Western blotting 1:1000-2000 ICH 1:1000-2000


Expiration Date

12 months from date of receipt.



Microtiter wells were coated with human ...

Sandwich ELISA



A sandwich ELISA was performed using ant...



Western blot analysis of extracts from H...