



## Product Datasheet

### SARS-CoV-2 Nucleoprotein-20 uL (orb1945984)

|                            |                                                                                                                 |
|----------------------------|-----------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | SARS-CoV-2 Nucleoprotein-20 uL                                                                                  |
| <b>Reactivity</b>          | Other                                                                                                           |
| <b>Tested Applications</b> | ELISA                                                                                                           |
| <b>Preservatives</b>       | 1X PBS, 0.02% NaN <sub>3</sub>                                                                                  |
| <b>Concentration</b>       | 0.5 mg/ml                                                                                                       |
| <b>Storage</b>             | 2-8°C                                                                                                           |
| <b>Tag</b>                 | His tag                                                                                                         |
| <b>Note</b>                | For research use only                                                                                           |
| <b>Isotype</b>             | SARS-CoV-2                                                                                                      |
| <b>Uniprot ID</b>          | <b>P0DTC9</b>                                                                                                   |
| <b>Dilution Range</b>      | Control antigen for SARS-CoV-2 (COVID-19) diagnostic assays. Contains C-terminal His tag used for purification. |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                  |



SDS-PAGE gel shows high purity for SARS-...



Microtiter wells were coated with SARS-C...