



Product Datasheet

SARS-CoV-2 NP (84D7) rabbit mAb Biotin Conjugate Antibody (orb1945987)

Description	SARS-CoV-2 NP (84D7) rabbit mAb Biotin Conjugate Antibody
Reactivity	Other
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	SARS-CoV-2 nucleoprotein (NP) specific peptide
Preservatives	1X PBS, 0.02% NaN ₃
Concentration	0.5 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	P0DTC9
Dilution Range	Used for specific, high-sensitivity detection of SARS-CoV-2 nucleocapsid protein (NP) in immunoassay. Can be paired with other NP-specific clones for detection in sandwich ELISA format.
Expiration Date	12 months from date of receipt.



Microtiter wells were coated with SARS-CoV-2 NP (84D7) rabbit mAb Biotin Conjugate Antibody.



A sandwich ELISA was performed using SARS-CoV-2 NP (84D7) rabbit mAb Biotin Conjugate Antibody.