

Product Datasheet

**Phospho-PKC (pan) (gamma Thr514) (PF4) rabbit
mAb Antibody (orb1946008)**

Description

Phospho-PKC (pan) (gamma Thr514) (PF4)

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC, WB

Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Thr514 of human phospho PKC γ

Preservatives

1X PBS, 0.02% NaN₃, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality

Monoclonal

Uniprot ID


P05771, **Q05655**, **P17252**, **P24723**, **Q04759**, **P05129**, **Q02156**

Dilution Range


1 μ g/mL – 0.001 μ g/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Expiration Date


12 months from date of receipt.



Flow cytometric analysis of HT1080 cells...



Flow cytometric analysis of Jurkat cells...



Peptide blocking flow cytometric analysis...