

Product Datasheet

Phospho-MCM2 (Ser139) (B12) rabbit mAb Antibody (orb1946023)

Description

Phospho-MCM2 (Ser139) (B12) rabbit mAb

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Ser139 of human phospho MCM2

Preservatives

1X PBS, 0.02% NaN₃, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality

Monoclonal

Uniprot ID


P49736

Dilution Range


1µg/mL – 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Expiration Date


12 months from date of receipt.



Flow
cytometric
analysis of
C6 cells,
se...



Western
blot
analysis of
L929 cell
extra...



Peptide
blocking
flow
cytometric
analysi...