

Product Datasheet

Phospho-Lyn (Tyr507) (5B6) rabbit mAb Antibody (orb1946031)

Description

Phospho-Lyn (Tyr507) (5B6) rabbit mAb Antibody

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

Applications

FC, WB

Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Tyr507 of human phospho Lyn

Preservatives

1X PBS, 0.02% NaN₃, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality

Monoclonal

Uniprot ID


P07948

Dilution Range


1µg/mL – 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Expiration Date


12 months from date of receipt.



Flow
cytometric
analysis of
Jurkat
cells...



Peptide
blocking
flow
cytometric
analysis...



Flow
cytometric
analysis of
C2C12
cells ...