


Product Datasheet

Phospho-PKCa (Thr497) (F1) rabbit mAb Antibody (orb1946091)


Description

Phospho-PKCa (Thr497) (F1) rabbit mAb


Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Thr497 of human phospho PKCα
Preservatives	1X PBS, 0.02% NaN ₃ , 50% Glycerol, 0.1% rAlbumin
Concentration	0.5 mg/ml
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgGκ
Clonality	Monoclonal
Uniprot ID	P17252
Dilution Range	1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
K562 cells
S...



Western
blot
analysis of
HT1080
cell ext...



Western
blot
analysis of
HT1080
cell ext...