

Product Datasheet

Cleaved PARP (Asp214) (H8) rabbit mAb Antibody (orb1946095)

Description

Cleaved PARP (Asp214) (H8) rabbit mAb

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, WB

Immunogen

A synthetic peptide corresponding to residues surrounding Asp214 of human PARP

Preservatives

1X PBS, 0.02% NaN₃, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality

Monoclonal

Uniprot ID


P09874

Dilution Range


1µg/mL – 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Expiration Date


12 months from date of receipt.



Western blot analysis of HeLa cell extra...



Peptide blocking flow cytometric analysis...



Flow cytometric analysis of SK.N.MC cell...