


Product Datasheet

Phospho-MEK1/2 (Ser221) (D3) rabbit mAb Antibody (orb1946141)


Description

Phospho-MEK1/2 (Ser221) (D3) rabbit mAb


Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Ser221 of human phospho MEK1/2.
Preservatives	1X PBS, 0.02% NaN ₃ , 50% Glycerol, 0.1% rAlbumin
Concentration	0.5 mg/ml
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	P36507 , Q02750
Dilution Range	1µg/mL – 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
Expiration Date	12 months from date of receipt.



Western blot analysis of HT1080 cell ext...



Western blot analysis of C2C12 cell extr...



Western blot analysis of NIH3T3 cell ext...