

Product Datasheet

Phospho-Histone H2A.X (Ser139) (1E4) rabbit mAb Antibody (orb1946178)

Description

Phospho-Histone H2A.X (Ser139) (1E4) rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

FC, WB

Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Ser139 of human phospho histone H2A.X.

Preservatives

1X PBS, 0.02% NaN₃, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGκ

Clonality

Monoclonal

Uniprot ID


P16104

Dilution Range


1µg/mL – 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Expiration Date


12 months from date of receipt.



Flow
cytometric
analysis of
3T3 cells
se...



Peptide
blocking
flow
cytometric
analysis...



Western
blot
analysis of
HEK293T
cell ex...