

## Product Datasheet

### Phospho-RelB (Ser552) (A7) rabbit mAb Antibody (orb1946198)

## Description

Phospho-RelB (Ser552) (A7) rabbit mAb Antibody

### Reactivity

Human, Mouse

### Conjugation

Unconjugated

### Tested

FC

### Applications

### Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Ser552 of human phospho RelB

### Preservatives

1X PBS, 0.02% NaN<sub>3</sub>, 50% Glycerol, 0.1% rAlbumin

### Concentration

0.5 mg/ml

### Storage

-20°C

### Note

For research use only

### Isotype

Rabbit IgGk

### Clonality

Monoclonal

### Uniprot ID


[Q01201](#)

### Dilution Range


1µg/mL – 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

### Expiration Date


12 months from date of receipt.



Flow  
cytometric  
analysis of  
NIH3T3  
cells...



Peptide  
blocking  
flow  
cytometric  
analysis...



Flow  
cytometric  
analysis of  
HeLa cells  
u...