

Product Datasheet

**Phospho-Chk1 (Ser345) (R3F9) rabbit mAb
Antibody (orb1946212)**

Description

Phospho-Chk1 (Ser345) (R3F9) rabbit mAb

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC

Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Ser345 of human phospho Chk1

Preservatives

1X PBS, 0.02% NaN₃, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality

Monoclonal

Uniprot ID


[O14757](#)

Dilution Range


1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Expiration Date


12 months from date of receipt.



Flow
cytometric
analysis of
Hela cells
S...



Peptide
blocking
flow
cytometric
analysis...



Chk1S345-
R3F9
recognizes
basal
phosphory...