


Product Datasheet

Phospho-MSK1 (Thr581) (A5) rabbit mAb Antibody (orb1946223)


Description

Phospho-MSK1 (Thr581) (A5) rabbit mAb


Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Thr581 of human phospho MSK1
Preservatives	1X PBS, 0.02% NaN ₃ , 50% Glycerol, 0.1% rAlbumin
Concentration	0.5 mg/ml
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	075582
Dilution Range	1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T cell extra...



Flow cytometric analysis of U937 cells S...



Peptide blocking flow cytometric analysis...