


Product Datasheet

Phospho-MKK7 (Ser271/Thr275) (R4F9) rabbit mAb APC conjugate Antibody (orb1946232)

Description	Phospho-MKK7 (Ser271/Thr275) (R4F9) rabbit mAb APC conjugate Antibody
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Ser271 and Thr275 of human MKK7
Preservatives	1X PBS, 0.09% NaN ₃ , 0.2% rAlbumin
Storage	2-8°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	O14733
Dilution Range	For flow cytometric staining, the suggested use of this reagent is 5 µL per million cells or 5 µL per 100 µL of staining volume. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
Expiration Date	12 months from date of receipt.



Flow
Cytometry
analysis
of 293T
cells un...