



Product Datasheet

Phospho-Chk2 (Thr68) (D12) rabbit mAb Antibody (orb1946260)


Description	Phospho-Chk2 (Thr68) (D12) rabbit mAb Antibody
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Thr68 of human phospho Chk2
Preservatives	1X PBS, 0.02% NaN ₃ , 50% Glycerol, 0.1% rAlbumin
Concentration	0.5 mg/ml
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	O96017
Dilution Range	1µg/mL – 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
HeLa cells t...



Peptide
blocking flow
cytometric
analysis...



Chk2T68-D12
recognizes
basal
phosphoryla...