

Product Datasheet

Phospho-CrkL (Tyr207) (G4) rabbit mAb Antibody (orb1946278)

Description

Phospho-CrkL (Tyr207) (G4) rabbit mAb Antibody

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Tyr207 of human phospho CrkL

Preservatives

1X PBS, 0.02% NaN₃, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality

Monoclonal

Uniprot ID


P46109

Dilution Range

1µg/mL – 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.


Expiration Date

12 months from date of receipt.




Phospho-CrkL (Y207)

Flow
cytometric
analysis of
K562 cells
t...



Western
blot
analysis of
K562 cell
extra...



Sample	Median	SD
Pre-CrkL	1712	1712
Post-CrkL	1712	1712
Pre-CrkL	1712	1712
Post-CrkL	1712	1712
Pre-CrkL	1712	1712
Post-CrkL	1712	1712
Pre-CrkL	1712	1712
Post-CrkL	1712	1712

Peptide
blocking
flow
cytometric
analysis...