

Product Datasheet

Phospho-ATF2 (Thr71) (G3) rabbit mAb Antibody (orb1946307)

Description

Phospho-ATF2 (Thr71) (G3) rabbit mAb Antibody

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

A synthetic phospho-peptide corresponding to residues surrounding Thr71 of human phospho ATF2

Preservatives

1X PBS, 0.02% NaN₃, 50% Glycerol, 0.1% rAlbumin

Concentration

0.5 mg/ml

Storage

-20°C

Note

For research use only

Isotype

Rabbit IgGk

Clonality

Monoclonal

Uniprot ID


P15336

Dilution Range


1µg/mL – 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Expiration Date


12 months from date of receipt.



Flow
cytometric
analysis of
NIH3T3
cells...



Western
blot
analysis of
293T cell
extra...



Peptide
blocking
flow
cytometric
analysis...