


Product Datasheet

Phospho-Histone H3 (Ser28) (D6) rabbit mAb Antibody (orb1946316)


Description

Phospho-Histone H3 (Ser28) (D6) rabbit mAb


Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Ser28 of human phospho Histone H3
Preservatives	1X PBS, 0.02% NaN ₃ , 50% Glycerol, 0.1% rAlbumin
Concentration	0.5 mg/ml
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	Q16695
Dilution Range	1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
NIH3T3
cells...



Western
blot
analysis of
HeLa cell
extra...



Flow
cytometric
analysis of
HeLa
seconda...