


Product Datasheet

Phospho-4E-BP1 (Thr37/46) (A5) rabbit mAb SureLight 488 conjugate Antibody (orb1946317)

Description	Phospho-4E-BP1 (Thr37/46) (A5) rabbit mAb SureLight 488 conjugate Antibody
Reactivity	Human, Mouse, Rat
Conjugation	SureLight 488
Tested Applications	FC
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Thr37/46 of human phospho 4E-BP1.
Preservatives	1X PBS, 0.09% NaN ₃ , 0.2% rAlbumin
Storage	2-8°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	Q13541
Dilution Range	For flow cytometric staining, the suggested use of this reagent is 5 µL per million cells or 5 µL per 100 µL of staining volume. It is recommended that the reagent be titrated for optimal performance for each application.
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Jurkat
cells...