



Product Datasheet

Phospho-4E-BP1 (Thr37/46) (A5) rabbit mAb Antibody (orb1946321)


Description	Phospho-4E-BP1 (Thr37/46) (A5) rabbit mAb Antibody
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Thr37/46 of human phospho 4E-BP1.
Preservatives	1X PBS, 0.02% NaN ₃ , 50% Glycerol, 0.1% rAlbumin
Concentration	0.5 mg/ml
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	Q13541
Dilution Range	1µg/mL – 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of Jurkat cells...



Immunohistochemical staining of rat kidn...



Western blot analysis of K562 cell extra...