


Product Datasheet

**Phospho-S6 Ribosomal Protein (Ser235/236)
(R3A2) rabbit mAb Antibody (orb1946327)**


Description

Phospho-S6 Ribosomal Protein (Ser235/236)


Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	A synthetic phospho-peptide corresponding to residues surrounding Ser235/236 of human phospho S6 Ribosomal Protein
Preservatives	1X PBS, 0.02% NaN3, 50% Glycerol, 0.1% rAlbumin
Concentration	0.5 mg/ml
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgGk
Clonality	Monoclonal
Uniprot ID	P62753
Dilution Range	1µg/mL - 0.001µg/mL. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH3T3 cell ext...



Flow cytometric analysis of U937 cells S...



Flow cytometric analysis of U937 cells, ...