

Product Datasheet

CHD5 antibody (orb19466)

Description

Goat polyclonal antibody to CHD5

Species/Host

Goat

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, IF

Applications
Target

CHD5

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:64000. WB: Preliminary experiments gave an approx 90kDa band in Human Brain (Cerebellum and Frontal Cortex) and Adrenal Gland lysates after 1 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 223kDa according to NP_056372.1. The 90kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

QSRRQLKSDRDKP

Clonality

Polyclonal

MW


223

Entrez
26038
Dilution Range


ELISA: 1:64000, IF/ICC: 10µg/ml, FACS: 10ug/ml

Expiration Date

12 months from date of receipt.



Immunofluorescence analysis of paraforma...



Flow cytometric analysis of paraformalde...