

Product Datasheet

CTDSP1 antibody (orb19467)

Description

Goat polyclonal antibody to CTDSP1

Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	CTDSP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Preliminary experiments gave an approx 60kDa band in human spleen, lung and placenta lysates after 0.3 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 29.2kDa according to NP_067021.1 and 29.1kDa according to NP_872580.1. The 60kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported
Protein Sequence	RGKGDQKSAASQK
Clonality	Polyclonal
MW	29.2; 29.1
Entrez	58190
Dilution Range	ELISA: 1:32000, WB: 0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.