


Product Datasheet

M6prbp1 antibody (orb19472)

Description	Goat polyclonal to M6prbp1
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	M6prbp1 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 48kDa band observed in Rat Skin lysates (calculated MW of 47.3kDa according to Mouse NP_080112.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	NQKTPQDAEKDPAK
Clonality	Polyclonal
MW	47.3
Dilution Range	ELISA: 1:32000, IF/ICC: 10µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH3T3 ce...