

## Product Datasheet

### Human HSD17 $\beta$ 10(17-Beta-Hydroxysteroid Dehydrogenase Type 10) ELISA Kit (orb1947263)

**Description**

Human HSD17β10(17-Beta-

10) ELISA Kit

**Reactivity**

Human

**Range**

25-1600pg/mL

**Target**


HSD17β10

**Note**

For research use only

**Application notes**

This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Human HSD17β10. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Human HSD17β10 and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Human HSD17β10, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm ± 2 nm. The OD value is proportional to the concentration of Human HSD17β10. You can calculate the concentration of Human HSD17β10 in the samples by comparing the OD of the samples to the standard curve.


**Sample Types**

Tissue homogenate

**Uniprot ID**
**Q99714**
**Sensitivity**

15 pg/mL

**Expiration Date**

Please enquire.