

Product Datasheet

FOXO4 antibody (orb19474)

Description	Goat polyclonal antibody to FOXO4
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FOXO4 / MLLT7
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 60-65kDa band observed in lysates of cell line K562 (calculated MW of 53.7kDa according to NP_005929.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	EAASQDRMPQDL
Clonality	Polyclonal
MW	53.7
Entrez	4303
Dilution Range	ELISA: 1:64000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 cell lysat...