



Product Datasheet

PSME2 antibody (orb19476)

Description	Goat polyclonal antibody to PSME2
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	PSME2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 27kDa band observed in Human Spleen lysates (calculated MW of 27.4kDa according to NP_002809.2). Recommended concentration: 0.5-1 µg/ml.
Protein Sequence	NLEKIVNPKGEEKP
Clonality	Polyclonal
MW	27.4
Entrez	5721
Dilution Range	ELISA: 1:32000, WB: 0.5-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb19476
(0.5ug/ml)
staining of
Human Sp...



orb19476
(5µg/ml)
staining of
paraffin ...