

Product Datasheet

PLK1 antibody (orb19485)

Description

PLK1 antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human PLK1. Please contact us for the exact immunogen sequence. The peptide is available as orb393969.

Target

PLK1

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW



68 kDa

Uniprot ID
P10683
NCBI
015973, 057057
Entrez
51083
Dilution Range


IHC-P: 1:100-600, WB: 1:200-1000

Expiration Date


12 months from date of receipt.



IHC-P image of human ovarian cancer tiss...



IHC-P staining of human breast tumour ti...



Immunohistochemical staining of human ov...