



Product Datasheet

SLC7A11 antibody (orb19488)


Description	Goat polyclonal antibody to SLC7A11
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	SLC7A11
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 50kDa band observed in Human Spleen and Human Thymus lysates (calculated MW of 55.4kDa according to NP_055146.1). Recommended concentration: 0.02-0.06 µg/ml.
Protein Sequence	KGQTQNFKDAFSGRD
Clonality	Polyclonal
MW	55.4
Entrez	23657
Dilution Range	ELISA: 1:32000, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cell...



Immunofluorescence analysis of A549 cell...



Immunohistochemical staining of Human Ki...