


Product Datasheet

Human GPR37(G Protein Coupled Receptor 37) ELISA Kit (orb1949074)

Description

Human GPR37(G Protein Coupled Receptor



Reactivity	Human
Range	0.16-10ng/mL
Target	GPR37
Note	For research use only
Application notes	This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Human GPR37. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Human GPR37 and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Human GPR37, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm $\pm$ 2 nm. The OD value is proportional to the concentration of Human GPR37. You can calculate the concentration of Human GPR37 in the samples by comparing the OD of the samples to the standard curve.
Sample Types	serum, plasma, Tissue homogenate and Other biological samples
Uniprot ID	O15354
Sensitivity	0.1 ng/mL
Expiration Date	Please enquire.