


Product Datasheet

arg1(rat) antibody (orb19491)

Description	Goat polyclonal to arg1(rat)
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Arginase, type 1 / arg1(rat)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide with sequence C-NHKPETDYLKPPK, from the C Terminus of the protein.ELISA: antibody detection limit dilution 1:128000.WB: Approx 37kDa band observed in Mouse and Rat Liver lysates (calculated MW of 34.8kDa according to mouse NP_031508.1). Recommended concentration: 0.03-0.1µg/ml.
Protein Sequence	NHKPETDYLKPPK
Clonality	Polyclonal
MW	35
Dilution Range	ELISA: 1:32000, WB: 0.01-0.03 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Liver lysat...



Western blot analysis of Mouse Liver (La...