

Product Datasheet

NR0B1 antibody (orb19492)

Description	Goat polyclonal antibody to NR0B1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	DAX1 / NR0B1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 50-55kDa band observed in human ovary and human testis lysates (calculated MW of 51.7kDa according to NP_000466.2). Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	CGEDHPQQGSTLY
Clonality	Polyclonal
MW	51.7
Entrez	190
Dilution Range	ELISA: 1:32000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



orb19492
(0.5ug/ml)
staining of
human te...