

Product Datasheet

Mouse ANXA2(Annexin A2) ELISA Kit (orb1949750)

Description	Mouse ANXA2(Annexin A2) ELISA Kit
Reactivity	Mouse
Range	1.57-100ng/mL
Target	ANXA2
Note	For research use only
Application notes	<p>This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Mouse ANXA2. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Mouse ANXA2 and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Mouse ANXA2, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm \pm 2 nm. The OD value is proportional to the concentration of Mouse ANXA2. You can calculate the concentration of Mouse ANXA2 in the samples by comparing the OD of the samples to the standard curve.</p>
Sample Types	Tissue homogenate and Cell lysate
Uniprot ID	P07356
Sensitivity	0.94 ng/mL
Expiration Date	Please enquire.

