


Product Datasheet

Goat Anti-TACR1 antibody (orb19502)

Description	Goat polyclonal antibody to Goat Anti-TACR1
Species/Host	Goat
Reactivity	Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Target	TACR1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 36kDa band observed in Human Testis lysates (calculated MW of 35.7kDa according to NP_056542.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	EIPGDSSDRYHEQ
Clonality	Polyclonal
MW	46.3; 35.7
Entrez	6869
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19502
(1µg/ml)
staining of
Human
Test...