

Product Datasheet

MDH2 Rabbit mAb Antibody (orb1950323)

Description	MDH2 Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P40926
Expiration Date	12 months from date of receipt.