


Product Datasheet

STRA6 antibody (orb19504)

Description	Goat polyclonal antibody to STRA6
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	STRA6 (mouse, rat)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 85kDa band observed in Rat Brain lysates (calculated MW of 73.8kDa according to NP_001025095.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	QAFRKAALTSAKAN
Clonality	Polyclonal
MW	73.7, 73.8
Dilution Range	ELISA: 1:32000, WB: 1-2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Brain lysat...



Western blot analysis of T3-L1 lysates u...